



SANITIZING BUCKET & SPRAY BOTTLE



SPRAY BOTTLE: 1/8 tsp. bleach
+ 2 cups water (75°F)
= 50 ppm to 100 ppm chlorine

(Label the bottle as to its contents)

BUCKET: 3/4 tsp. bleach
+ 1 gallon water (75°F)
= 50 ppm to 100 ppm chlorine

(Label the bucket as to its contents)

- **CHANGE THE SOLUTION FREQUENTLY!**
- **Use test strips to check the concentration**
- **Store the wiping cloths in the solution**

(Effectiveness is hindered by organic matter,
an alkaline (basic) pH and by high water temperature.)